

For Research Use Only

PROCA1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83193-4-PBS

Basic Information

Catalog Number: 83193-4-PBS	GenBank Accession Number: BC029574	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 147011	CloneNo.: 230538G2
Source: Rabbit	UNIPROT ID: Q8NCQ7	
Isotype: IgG	Full Name: proline-rich cyclin A1-interacting protein	
Immunogen Catalog Number: AG10851	Calculated MW: 336 aa, 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00294-2, PROCA1 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83193-3-PBS. Detection antibody: 83193-4-PBS. Standard: Ag10851. Range: 5-40 ng/mL.