

For Research Use Only

CEACAM6 Recombinant antibody, PBS Only



Catalog Number: 83203-5-PBS

Basic Information

Catalog Number:

83203-5-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005008

GeneID (NCBI):

4680

UNIPROT ID:

P40199

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 6 (non-
specific cross reacting antigen)

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230475D5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83203-5-RR (CEACAM6 antibody) at dilution of 1:7000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83203-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83203-5-RR against Human CEACAM6 were performed. The affinity constant is 0.113 nM.