

For Research Use Only

CEP350-Specific Recombinant antibody



Catalog Number: 83208-2-RR

Basic Information

Catalog Number: 83208-2-RR	GenBank Accession Number: NM_014810	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 9857	CloneNo.: 230320D4
Source: Rabbit	UNIPROT ID: Q5VT06	
Isotype: IgG	Full Name: centrosomal protein 350kDa	
	Calculated MW: 351 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

CEP350, also named as CAP350, FLJ38282, FLJ44058, GM133 and KIAA0480, plays an essential role in centriole growth by stabilizing a procentriolar seed composed of at least, SASS6 and CENPJ. It is required for anchoring microtubules to the centrosomes and for the integrity of the microtubule network. CEP350 recruits PPARA to discrete subcellular compartments and thereby modulates PPARA activity.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

