

For Research Use Only

SP8 Recombinant antibody

Catalog Number: 83230-3-RR



Basic Information

Catalog Number: 83230-3-RR	GenBank Accession Number: BC038669	Purification Method: Protein A purification
Size: 100ul , Concentration: 300 µg/ml by Nanodrop;	GeneID (NCBI): 221833	CloneNo.: 240011G4
Source: Rabbit	UNIPROT ID: Q8IXZ3	
Isotype: IgG	Full Name: Sp8 transcription factor	
Immunogen Catalog Number: AG22934	Calculated MW: 508 aa, 51 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83230-3-RR against Human SP8 were performed. The affinity constant is 5.44 nM.