

For Research Use Only

# DUS2L Recombinant antibody, PBS Only (Capture)

Catalog Number: 83265-1-PBS



## Basic Information

Catalog Number:

83265-1-PBS

Size:

100ug , Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23831

GenBank Accession Number:

BC006527

GeneID (NCBI):

54920

UNIPROT ID:

Q9NX74

Full Name:

dihydrouridine synthase 2-like, SMM1  
homolog (S. cerevisiae)

Purification Method:

Protein A purification

CloneNo.:

240192A6

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

# Selected Validation Data



Standard curve of MP00261-1, Human DUS2L  
Recombinant Matched Antibody Pair - BSA and  
Azide Free measured by Cytometric bead array.  
Capture antibody: 83265-1-PBS. Detection  
antibody: 83265-4-PBS. Standard: Ag23831. Range:  
1.25-40 ng/mL