

For Research Use Only

DUS2L Recombinant antibody, PBS Only (Detector)

Catalog Number: 83265-4-PBS



Basic Information

Catalog Number: 83265-4-PBS	GenBank Accession Number: BC006527	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 54920	CloneNo.: 240192B1
Source: Rabbit	UNIPROT ID: Q9NX74	
Isotype: IgG	Full Name: dihydrouridine synthase 2-like, SMM1 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG23831		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83265-4-PBS targets DUS2L as part of a matched antibody pair.

MP00261-1: 83265-3-PBS capture and 83265-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00261-1, DUS2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83265-3-PBS. Detection antibody: 83265-4-PBS. Standard: Ag23831. Range: 1.25-40 ng/mL.