

For Research Use Only

Mouse SHH Recombinant antibody, PBS Only



Catalog Number: 83280-1-PBS

Basic Information

Catalog Number: 83280-1-PBS	GenBank Accession Number: NM_009170.3	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 20423	CloneNo.: 240057E7
Source: Rabbit	UNIPROT ID: Q62226	
Isotype: IgG	Full Name: sonic hedgehog	
	Calculated MW: 48KD	

Applications

Tested Applications:
WB, Indirect ELISA
Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83280-1-RR (SHH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83280-1-PBS in a different storage buffer formulation.