

For Research Use Only

Mouse Albumin Recombinant antibody, PBS Only (Detector)



Catalog Number: 83289-2-PBS

Basic Information

Catalog Number: 83289-2-PBS	GenBank Accession Number: NM_009654.4	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 11657	CloneNo.: 240245H4
Source: Rabbit	UNIPROT ID: P07724	
Isotype: IgG	Full Name: albumin	
	Calculated MW: 69 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mouse

Product Information

83289-2-PBS targets Albumin as part of a matched antibody pair:

MP00329-2: 83289-1-PBS capture and 83289-2-PBS detection (validated in Cytometric bead array)

MP00329-3: 83289-3-PBS capture and 83289-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00329-2, MOUSE Albumin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83289-1-PBS. Detection antibody: 83289-2-PBS. Standard: Eg0886. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00329-3, MOUSE Albumin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83289-3-PBS. Detection antibody: 83289-2-PBS. Standard: Eg0886. Range: 0.313-40 ng/mL.