

For Research Use Only

# Mouse Albumin Recombinant antibody, PBS Only



Catalog Number: 83289-4-PBS

## Basic Information

Catalog Number:

83289-4-PBS

Size:

100ug, Concentration: 1 mg/ml by  
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_009654.4

GeneID (NCBI):

11657

UNIPROT ID:

P07724

Full Name:

albumin

Calculated MW:

69 kDa

Purification Method:

Protein A purification

CloneNo.:

240245E1

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Mouse

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83289-4-RR (Albumin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83289-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83289-4-RR against Mouse Albumin were performed. The affinity constant is 0.883 nM.