

For Research Use Only

CRIP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83340-4-PBS



Basic Information

| | | |
|---|---|---|
| Catalog Number: 83340-4-PBS | GenBank Accession Number: BC002738 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 1396 | CloneNo.: 240176D8 |
| Source: Rabbit | UNIPROT ID: P50238 | |
| Isotype: IgG | Full Name: cysteine-rich protein 1 (intestinal) | |
| Immunogen Catalog Number: AG7588 | Calculated MW: 9 kDa | |

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

83340-4-PBS targets CRIP1 as part of a matched antibody pair:

MP00382-3: 83340-4-PBS capture and 83340-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00382-3, Human CRIP1 Recombinant Matched Antibody Pair - PBS only. 83340-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7588. 83340-3-PBS was HRP conjugated as the detection antibody. Range: 46.88-3000 pg/mL.



The mean CRIP1 concentration was determined to be 3,364.0 pg/mL in PC-3 cell extract based on a 3.0 mg/mL extract load and 3,015.0 pg/mL in HeLa cell extract based on a 3.0 mg/mL extract load.