

MTF1 Recombinant antibody, PBS Only (Capture)

Catalog Number: **83344-2-PBS**

Basic Information

| | | |
|--|--|---|
| Catalog Number: 83344-2-PBS | GenBank Accession Number: BC014454 | Purification Method: Protein A purification |
| Size: 100ug, Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 4520 | CloneNo.: 240336H1 |
| Source: Rabbit | UNIPROT ID: Q14872 | |
| Isotype: IgG | Full Name: metal-regulatory transcription factor 1 | |
| Immunogen Catalog Number: AG21848 | Calculated MW: 753 aa, 81 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83344-2-PBS targets MTF1 as part of a matched antibody pair:

MP00386-2: 83344-2-PBS capture and 83344-3-PBS detection (validated in Cytometric bead array)

MP00386-3: 83344-2-PBS capture and 83344-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00386-2, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00386-3, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-2-PBS. Detection antibody: 83344-1-PBS. Standard: Ag21848. Range: 0.625-20 ng/mL