

For Research Use Only

BAP1 Recombinant antibody

Catalog Number: 83371-1-RR



Basic Information

Catalog Number:

83371-1-RR

Size:

100ul , Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34415

GenBank Accession Number:

BC001596

GeneID (NCBI):

8314

UNIPROT ID:

Q92560

Full Name:

BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase)

Calculated MW:

81-91 kDa

Observed MW:

91-100 kDa

Purification Method:

Protein A purification

CloneNo.:

240349A7

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HEK-293 cells, LNCaP cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83371-1-RR (BAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83371-1-RR against Human BAP1 were performed. The affinity constant is below 1 pM.