

For Research Use Only

DSC3 Recombinant antibody, PBS Only



Catalog Number: 83386-5-PBS

Basic Information

Catalog Number:

83386-5-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34648

GenBank Accession Number:

NM_001941

GeneID (NCBI):

1825

UNIPROT ID:

Q14574

Full Name:

desmocollin 3

Calculated MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

240292B12

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com




W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83386-5-RR (DSC3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83386-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83386-5-RR against Human DSC3 were performed. The affinity constant is below 1 pM.