

For Research Use Only

# NOL11 Recombinant antibody

Catalog Number: 83391-7-RR



## Basic Information

Catalog Number:

83391-7-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11783

GenBank Accession Number:

BC064404

GeneID (NCBI):

25926

UNIPROT ID:

Q9H8H0

Full Name:

nucleolar protein 11

Calculated MW:

719 aa, 81 kDa

Observed MW:

81 kDa

Purification Method:

Protein A purification

CloneNo.:

240371B9

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

