

For Research Use Only

MOSPD2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83398-2-PBS

Basic Information

Catalog Number: 83398-2-PBS	GenBank Accession Number: BC030641	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 158747	CloneNo.: 240282H6
Source: Rabbit	UNIPROT ID: Q8NHP6	
Isotype: IgG	Full Name: motile sperm domain containing 2	
Immunogen Catalog Number: AG22206	Calculated MW: 518 aa, 60 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83398-2-PBS targets MOSPD2 as part of a matched antibody pair:

MP00419-1: 83398-3-PBS capture and 83398-2-PBS detection (validated in Cytometric bead array)

MP00419-2: 83398-1-PBS capture and 83398-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00419-1, MOSPD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83398-3-PBS. Detection antibody: 83398-2-PBS. Standard: Ag22206. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00419-2, MOSPD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83398-1-PBS. Detection antibody: 83398-2-PBS. Standard: Ag22206. Range: 0.78-100 ng/mL