

For Research Use Only

KRT77 Recombinant antibody, PBS Only

Catalog Number: 83415-3-PBS



Basic Information

Catalog Number:

83415-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21942

GenBank Accession Number:

BC122558

GeneID (NCBI):

374454

UNIPROT ID:

Q7Z794

Full Name:

keratin 77

Calculated MW:

578 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

240490A11

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 83415-3-RR (KRT77 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83415-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83415-3-RR against Human KRT77 were performed. The affinity constant is 0.209 nM.