

For Research Use Only

Mouse CD117/c-Kit Recombinant antibody, PBS Only (Detector)

Catalog Number: 83419-2-PBS



Basic Information

Catalog Number: 83419-2-PBS	GenBank Accession Number: NM_001122733.1	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 16590	CloneNo.: 240246G11
Source: Rabbit	UNIPROT ID: P05532-1	
Isotype: IgG	Full Name: kit oncogene	
Immunogen Catalog Number: EG0883	Calculated MW: 109 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

83419-2-PBS targets CD117/c-Kit as part of a matched antibody pair:

MP00418-1: 83419-1-PBS capture and 83419-2-PBS detection (validated in Cytometric bead array)

MP00418-3: 83419-3-PBS capture and 83419-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00418-1, MOUSE CD117/c-kit Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83419-1-PBS. Detection antibody: 83419-2-PBS. Standard: Eg0883. Range: 0.781-25 ng/mL



Cytometric bead array standard curve of MP00418-3, MOUSE CD117/c-kit Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83419-3-PBS. Detection antibody: 83419-2-PBS. Standard: Eg0883. Range: 6.25-100 ng/mL