

For Research Use Only

LBP Recombinant antibody, PBS Only

Catalog Number: 83427-1-PBS



Basic Information

Catalog Number:

83427-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004139.4

GeneID (NCBI):

3929

UNIPROT ID:

P18428

Full Name:

Lipopolysaccharide binding protein

Calculated MW:

53KD

Purification Method:

Protein A purification

CloneNo.:

230524G10

Applications

Tested Applications:

ELISA, WB

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 83427-1-RR (LBP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83427-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83427-1-RR against Human LBP were performed. The affinity constant is 27.1 nM.