

For Research Use Only

THOC5 Recombinant antibody, PBS Only

Catalog Number: 83435-1-PBS



Basic Information

Catalog Number:

83435-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6664

GenBank Accession Number:

BC003615

GeneID (NCBI):

8563

UNIPROT ID:

Q13769

Full Name:

THO complex 5

Calculated MW:

79 kDa

Purification Method:

Protein A purification

CloneNo.:

240013F7

Applications

Tested Applications:

IF/ICC, FC (Intra), ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ U2OS cells were intracellularly stained with 0.25 ug Thoc5 Recombinant Antibody (83435-1-RR, Clone:240013F7) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). This data was developed using the same antibody clone with 83435-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83435-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed U2OS cells using THOC5 antibody (83435-1-RR, Clone: 240013F7) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83435-1-PBS in a different storage buffer formulation.