

For Research Use Only

LPAR4 Recombinant antibody, PBS Only



Catalog Number: 83457-1-PBS

Basic Information

Catalog Number: 83457-1-PBS	GenBank Accession Number: BC069996	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2846	CloneNo.: 240079C3
Source: Rabbit	UNIPROT ID: Q99677	
Isotype: IgG	Full Name: lysophosphatidic acid receptor 4	
Immunogen Catalog Number: AG17817	Calculated MW: 370 aa, 42 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83457-1-RR against Human LPAR4 were performed. The affinity constant is below 1 pM.