

For Research Use Only

ENT2 Recombinant antibody

Catalog Number: 83461-1-RR



Basic Information

Catalog Number:

83461-1-RR

GenBank Accession Number:

BC093634

Purification Method:

Protein A purification

Size:

100ul, Concentration: 1000 µg/ml by Nanodrop;

GeneID (NCBI):

3177

CloneNo.:

240090H1

Source:

Rabbit

UNIPROT ID:

Q14542

Full Name:

solute carrier family 29 (nucleoside transporters), member 2

Isotype:

IgG

Immunogen Catalog Number:

AG23310

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83461-1-RR against Human SLC29A2/ENT2 were performed. The affinity constant is 1.39 nM.