

For Research Use Only

# GSTT2B Recombinant antibody, PBS Only



Catalog Number: 83469-4-PBS

## Basic Information

<b>Catalog Number:</b> 83469-4-PBS	<b>GenBank Accession Number:</b> BC071700	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 653689	<b>CloneNo.:</b> 240493G4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P0CG30	
<b>Isotype:</b> IgG	<b>Full Name:</b> glutathione S-transferase theta 2B (gene/pseudogene)	
<b>Immunogen Catalog Number:</b> AG11828	<b>Calculated MW:</b> 244 aa, 28 kDa	

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
Human

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83469-4-RR (GSTT2B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83469-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83469-4-RR against Human GSTT2B were performed. The affinity constant is below 1 pM.