

For Research Use Only

# IL-2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 83558-7-PBS



## Basic Information

<b>Catalog Number:</b> 83558-7-PBS	<b>GenBank Accession Number:</b> GeneID (NCBI): 3558	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>UNIPROT ID:</b> P60568	<b>CloneNo.:</b> 252448F9
<b>Source:</b> Rabbit	<b>Full Name:</b> interleukin 2	
<b>Isotype:</b> IgG		
<b>Immunogen Catalog Number:</b> HZ-1015		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

83558-7-PBS targets IL-2 as part of a matched antibody pair:

MP00545-4: 83558-7-PBS capture and 83558-6-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS only, pH7.3


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Human peripheral blood mononuclear cells (PBMC) ( $1 \times 10^6$  cells/mL) were cultured in RPMI-1640 supplemented with 10% fetal bovine serum, 100 U/mL penicillin and 100  $\mu$ g/mL streptomycin sulfate. The cell culture supernatants were stimulated or unstimulated with 10  $\mu$ g/mL of PHA for 5 days. Aliquots of the cell culture supernates were removed and assayed for levels of human IL-2, and measured 7,841.4 pg/mL, 44.0 pg/mL respectively.



Cytometric bead array standard curve of MP00545-4, IL-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83558-7-PBS. Detection antibody: 83558-6-PBS. Standard: HZ-1015. Range: 7.81-1000 pg/mL.