

For Research Use Only

# Rat VEGFA Recombinant antibody, PBS Only

Catalog Number: 83610-3-PBS



## Basic Information

Catalog Number:

83610-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0941

GenBank Accession Number:

NM\_001287108.1

GeneID (NCBI):

83785

UNIPROT ID:

P16612-2

Full Name:

vascular endothelial growth factor A

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

240573H8

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

rat

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

