

For Research Use Only

MTMR15 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83622-4-PBS



Basic Information

Catalog Number: 83622-4-PBS	GenBank Accession Number: BC047882	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 22909	CloneNo.: 240629D12
Source: Rabbit	UNIPROT ID: Q9Y2M0	
Isotype: IgG	Full Name: myotubularin related protein 15	
Immunogen Catalog Number: AG11651	Calculated MW: 60 kDa, 114 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83622-4-PBS targets MTMR15 as part of a matched antibody pair:

MP00572-1: 83622-4-PBS capture and 83622-1-PBS detection (validated in Cytometric bead array)

MP00572-3: 83622-2-PBS capture and 83622-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00572-1, MTMR15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83622-4-PBS. Detection antibody: 83622-1-PBS. Standard: Ag11651. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00572-3, MTMR15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83622-2-PBS. Detection antibody: 83622-4-PBS. Standard: Ag11651. Range: 0.78-100 ng/mL