

For Research Use Only

LPXN Recombinant antibody, PBS Only



Catalog Number: 83657-5-PBS

Basic Information

Catalog Number: 83657-5-PBS	GenBank Accession Number: BC019035	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9404	CloneNo.: 240483G10
Source: Rabbit	UNIPROT ID: O60711	
Isotype: IgG	Full Name: leupaxin	
Immunogen Catalog Number: AG34956	Calculated MW: 43 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83657-5-RR (LPXN antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83657-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83657-5-RR against Human LPXN were performed. The affinity constant is 0.187 nM.