

For Research Use Only

Mouse Il-7 Recombinant antibody, PBS Only



Catalog Number: 83672-1-PBS

Basic Information

Catalog Number:

83672-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008371.5

GeneID (NCBI):

16196

UNIPROT ID:

Q544C8

Full Name:

interleukin 7

Calculated MW:

18kd

Purification Method:

Protein A purification

CloneNo.:

240343C8

Applications

Tested Applications:

ELISA, WB

Species Specificity:

mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Mouse IL-7 Recombinant proteins (Eg0403) and human IL-7 Recombinant protein (Eg0830) were subjected to SDS PAGE followed by western blot with 83672-1-RR (IL-7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83672-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83672-1-RR against Mouse IL-7 were performed. The affinity constant is 3.70 nM.