

For Research Use Only



Interferon alpha 1/IFNA1

Recombinant antibody, PBS Only

Catalog Number: 83684-7-PBS

Basic Information

Catalog Number:

83684-7-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_024013.3

GeneID (NCBI):

3447

UNIPROT ID:

P01562

Full Name:

interferon, alpha 13

Calculated MW:

22kDa

Purification Method:

Protein A purification

CloneNo.:

240667F11

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed LPS and Brefeldin A treated RAW 264.7 cells using Interferon alpha 1/IFNA1 antibody (83684-7-RR, Clone: 240667F11) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83684-7-PBS in a different storage buffer formulation.