

For Research Use Only

RABGAP1L Recombinant antibody, PBS Only



Catalog Number: 83698-4-PBS

Basic Information

Catalog Number:

83698-4-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4856

GenBank Accession Number:

BC041888

GeneID (NCBI):

9910

UNIPROT ID:

B7ZAP0

Full Name:

RAB GTPase activating protein 1-like

Calculated MW:

815 aa, 93 kDa

Purification Method:

Protein A purification

CloneNo.:

240663G12

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83698-4-RR (RABGAP1L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83698-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83698-4-RR against Human RABGAP1L were performed. The affinity constant is below 1 pM.