

For Research Use Only

GBP4 Recombinant antibody, PBS Only



Catalog Number: 83699-4-PBS

Basic Information

Catalog Number: 83699-4-PBS	GenBank Accession Number: BC070055	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 115361	CloneNo.: 240730H5
Source: Rabbit	UNIPROT ID: Q96PP9	
Isotype: IgG	Full Name: guanylate binding protein 4	
Immunogen Catalog Number: AG12006	Calculated MW: 640 aa, 73 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 83699-4-RR (GBP4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83699-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83699-4-RR against Human GBP4 were performed. The affinity constant is 11.6 pM.