

For Research Use Only

Rat Mmp9 Recombinant antibody, PBS Only



Catalog Number: 83768-3-PBS

Basic Information

Catalog Number:

83768-3-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_031055.2

GeneID (NCBI):

81687

UNIPROT ID:

P50282

Full Name:

matrix metalloproteinase 9

Calculated MW:

79kd

Purification Method:

Protein A purification

CloneNo.:

240895G9

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83768-3-RR (Mmp9 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83768-3-PBS in a different storage buffer formulation.