

For Research Use Only

GATA3 Recombinant antibody, PBS Only



Catalog Number: 83808-2-PBS

Basic Information

Catalog Number: 83808-2-PBS	GenBank Accession Number: BC003070	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2625	CloneNo.: 24068488
Source: Rabbit	UNIPROT ID: P23771	
Isotype: IgG	Full Name: GATA binding protein 3	
Immunogen Catalog Number: AG0174	Calculated MW: 443 aa, 48 kDa	

Applications

Tested Applications:
ELISA, IF/ICC

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using GATA3 antibody (83808-2-RR, Clone: 240684B8) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83808-2-PBS in a different storage buffer formulation.