

For Research Use Only

PSENEN Recombinant antibody, PBS Only



Catalog Number: 83818-1-PBS

Basic Information

Catalog Number: 83818-1-PBS	GenBank Accession Number: BC009575	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55851	CloneNo.: 240596H2
Source: Rabbit	UNIPROT ID: Q9NZ42	
Isotype: IgG	Full Name: presenilin enhancer 2 homolog (C. elegans)	
Immunogen Catalog Number: AG9243	Calculated MW: 101 aa, 12 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83818-1-RR (PSENEN antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83818-1-PBS in a different storage buffer formulation.