

For Research Use Only

FAM127B Recombinant antibody, PBS Only



Catalog Number: 83862-2-PBS

Basic Information

Catalog Number:

83862-2-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13930

GenBank Accession Number:

BC000393

GeneID (NCBI):

26071

UNIPROT ID:

Q9BWD3

Full Name:

family with sequence similarity 127,
member B

Calculated MW:

113 aa, 13 kDa

Purification Method:

Protein A purification

CloneNo.:

240716B7

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using FAM127B antibody (83862-2-RR, Clone: 240716B7) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83862-2-PBS in a different storage buffer formulation.