

For Research Use Only

CAMSAP1 Recombinant antibody, PBS Only



Catalog Number: 83864-1-PBS

Basic Information

Catalog Number:

83864-1-PBS

Size:

100ug , Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC144082

GeneID (NCBI):

157922

UNIPROT ID:

Q5T5Y3

Full Name:

calmodulin regulated spectrin-
associated protein 1

Calculated MW:

179 kDa

Purification Method:

Protein A purification

CloneNo.:

240633H1

Applications

Tested Applications:

ELISA, WB

Species Specificity:

mouse, human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83864-1-RR (CAMSAP1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83864-1-PBS in a different storage buffer formulation.