

For Research Use Only

# CRNKL1 Recombinant antibody, PBS Only



Catalog Number: 83914-1-PBS

## Basic Information

Catalog Number:

83914-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_016652

GeneID (NCBI):

51340

UNIPROT ID:

Q9BZJO

Full Name:

crooked neck pre-mRNA splicing factor-like 1 (Drosophila)

Calculated MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

240697D12

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CRNKL1 antibody (83914-1-RR, Clone: 240697D12 ) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83914-1-PBS in a different storage buffer formulation.