

For Research Use Only

FAM185A Recombinant antibody, PBS Only (Capture)



Catalog Number: 83921-1-PBS

Basic Information

Catalog Number: 83921-1-PBS	GenBank Accession Number: BC029175	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 222234	CloneNo.: 240804D12
Source: Rabbit	UNIPROT ID: Q8NOU4	
Isotype: IgG	Full Name: family with sequence similarity 185, member A	
Immunogen Catalog Number: AG15041	Calculated MW: 393 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83921-1-PBS targets FAM185A as part of a matched antibody pair:

MP00824-1: 83921-1-PBS capture and 83921-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00824-1, FAM185A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83921-1-PBS. Detection antibody: 83921-2-PBS. Standard: Ag15041. Range: 3.125-100 ng/mL