

For Research Use Only

CXXC5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83940-2-PBS

Basic Information

Catalog Number: 83940-2-PBS	GenBank Accession Number: BC017439	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51523	CloneNo.: 24103487
Source: Rabbit	UNIPROT ID: Q7LFL8	
Isotype: IgG	Full Name: CXXC finger 5	
Immunogen Catalog Number: AG9733	Calculated MW: 322aa,33 kDa; 227aa,24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00925-2, CXXC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83940-4-PBS. Detection antibody: 83940-2-PBS. Standard: Ag9733. Range: 0.313-40 ng/mL