

TLE3 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83955-2-PBS

Basic Information

Catalog Number: 83955-2-PBS	GenBank Accession Number: BC015729	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7090	CloneNo.: 240680D11
Source: Rabbit	UNIPROT ID: Q04726	
Isotype: IgG	Full Name: transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)	
Immunogen Catalog Number: AG1930	Calculated MW: 772 aa, 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83955-2-PBS targets TLE3 as part of a matched antibody pair:

MP00909-2: 83955-2-PBS capture and 83955-1-PBS detection (validated in Cytometric bead array)

MP00909-3: 83955-3-PBS capture and 83955-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00909-2, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-2-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00909-3, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-3-PBS. Detection antibody: 83955-2-PBS. Standard: Ag1930. Range: 0.313-20 ng/mL