

For Research Use Only

AGXT2L1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83984-4-PBS

Basic Information

Catalog Number: 83984-4-PBS	GenBank Accession Number: BC022526	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 64850	CloneNo.: 241119G4
Source: Rabbit	UNIPROT ID: Q8TBG4	
Isotype: IgG	Full Name: alanine-glyoxylate aminotransferase 2-like 1	
Immunogen Catalog Number: AG10584	Calculated MW: 499 aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83984-4-PBS targets AGXT2L1 as part of a matched antibody pair:

MP00885-2: 83984-2-PBS capture and 83984-4-PBS detection (validated in Cytometric bead array)

MP00885-3: 83984-4-PBS capture and 83984-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00885-2, AGXT2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83984-2-PBS. Detection antibody: 83984-4-PBS. Standard: Ag10584. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00885-3, AGXT2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83984-4-PBS. Detection antibody: 83984-3-PBS. Standard: Ag10584. Range: 0.625-80 ng/mL.