

For Research Use Only

NMI Recombinant monoclonal antibody

Catalog Number: 83986-1-RR



Basic Information

Catalog Number: 83986-1-RR	GenBank Accession Number: BC001268	Purification Method: Protein A purification
Size: 100ul	GeneID (NCBI): 9111	CloneNo.: 241135A3
Source: Rabbit	UNIPROT ID: Q13287	
Isotype: IgG	Full Name: N-myc (and STAT) interactor	
Immunogen Catalog Number: AG34759	Calculated MW: 35 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data