

For Research Use Only

CALCOCO1 Recombinant antibody

Catalog Number: 84009-5-RR



Basic Information

Catalog Number:

84009-5-RR

Size:

100ul, Concentration: 1000 ug/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13876

GenBank Accession Number:

BC003177

GeneID (NCBI):

57658

UNIPROT ID:

Q9P1Z2

Full Name:

calcium binding and coiled-coil domain 1

Calculated MW:

691 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

241127B3

Recommended Dilutions:

IF/ICC 1:125-1:500

Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : HT-1080 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HT-1080 cells using CALCOCO1 antibody (84009-5-RR, Clone: 241127B3) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red).



Biolayer interferometry (BLI) kinetic assays of 84009-5-RR against Human CALCOCO1 were performed. The affinity constant is below 1 pM.