

For Research Use Only

Siglec-8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84026-2-PBS



Basic Information

Catalog Number: 84026-2-PBS	GenBank Accession Number: BC053319	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 27181	CloneNo.: 241276A8
Source: Rabbit	UNIPROT ID: Q9NYZ4	
Isotype: IgG	Full Name: sialic acid binding Ig-like lectin 8	
	Calculated MW: 499 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

84026-2-PBS targets Siglec-8 as part of a matched antibody pair:

MP00960-3: 84026-2-PBS capture and 84026-4-PBS detection (validated in Cytometric bead array)

MP00960-4: 84026-2-PBS capture and 84026-1-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00960-4, Human Siglec-8 Recombinant Matched Antibody Pair - PBS only. 84026-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard RP02311. 84026-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



Serum of six individual healthy human donors was measured. The Siglec-8 concentration of detected samples was determined to be 7.56 ng/mL with a range of 0.19-32.32 ng/mL



Cytometric bead array standard curve of MP00960-3, SIGLEC8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84026-2-PBS. Detection antibody: 84026-4-PBS. Standard: RP02311. Range: 0.313-40 ng/mL