

For Research Use Only

# Siglec-8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84026-4-PBS



## Basic Information

<b>Catalog Number:</b> 84026-4-PBS	<b>GenBank Accession Number:</b> BC053319	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 27181	<b>CloneNo.:</b> 241276C1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9NYZ4	
<b>Isotype:</b> IgG	<b>Full Name:</b> sialic acid binding Ig-like lectin 8	
	<b>Calculated MW:</b> 499 aa, 54 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

84026-4-PBS targets Siglec-8 as part of a matched antibody pair:

MP00960-1: 84026-3-PBS capture and 84026-4-PBS detection (validated in Cytometric bead array)

MP00960-2: 84026-1-PBS capture and 84026-4-PBS detection (validated in Cytometric bead array)

MP00960-3: 84026-2-PBS capture and 84026-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


# Selected Validation Data



Cytometric bead array standard curve of MP00960-1, SIGLEC8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84026-3-PBS. Detection antibody: 84026-4-PBS. Standard: RP02311. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00960-2, SIGLEC8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84026-1-PBS. Detection antibody: 84026-4-PBS. Standard: RP02311. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00960-3, SIGLEC8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84026-2-PBS. Detection antibody: 84026-4-PBS. Standard: RP02311. Range: 0.313-40 ng/mL