

For Research Use Only

# PEG3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84030-4-PBS

## Basic Information

<b>Catalog Number:</b> 84030-4-PBS	<b>GenBank Accession Number:</b> BC150272	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5178	<b>CloneNo.:</b> 241234G1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9GZU2	
<b>Isotype:</b> IgG	<b>Full Name:</b> paternally expressed 3	
<b>Immunogen Catalog Number:</b> AG20132	<b>Calculated MW:</b> 1588 aa, 181 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

84030-4-PBS targets PEG3 as part of a matched antibody pair:

MP00940-2: 84030-2-PBS capture and 84030-4-PBS detection (validated in Cytometric bead array)

MP00940-3: 84030-1-PBS capture and 84030-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Cytometric bead array standard curve of MP00940-2, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-2-PBS. Detection antibody: 84030-4-PBS. Standard:null. Range: null.



Cytometric bead array standard curve of MP00940-2, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-2-PBS. Detection antibody: 84030-4-PBS. Standard: Ag20132. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP00940-3, PEG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84030-1-PBS. Detection antibody: 84030-4-PBS. Standard: Ag20132. Range: 0.156-20 ng/mL.