

For Research Use Only

PCF11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84081-3-PBS



Basic Information

Catalog Number: 84081-3-PBS	GenBank Accession Number: BC146778	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51585	CloneNo.: 241213E8
Source: Rabbit	UNIPROT ID: O94913	
Isotype: IgG	Full Name: PCF11, cleavage and polyadenylation factor subunit, homolog (S. cerevisiae)	
Immunogen Catalog Number: AG20263	Calculated MW: 1555 aa, 173 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84081-3-PBS targets PCF11 as part of a matched antibody pair:

MP00996-1: 84081-2-PBS capture and 84081-3-PBS detection (validated in Cytometric bead array)

MP00996-2: 84081-1-PBS capture and 84081-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00996-1, PCF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84081-2-PBS. Detection antibody: 84081-3-PBS. Standard: Ag20263. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00996-2, PCF11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84081-1-PBS. Detection antibody: 84081-3-PBS. Standard: Ag20263. Range: 0.625-80 ng/mL