

For Research Use Only

CD244 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84086-2-PBS



Basic Information

Catalog Number: 84086-2-PBS	GenBank Accession Number: BC053985	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 51744	CloneNo.: 241260E3
Source: Rabbit	UNIPROT ID: Q9BZW8	
Isotype: IgG	Full Name: CD244 molecule, natural killer cell receptor 2B4	
	Calculated MW: 365 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84086-2-PBS targets CD244 as part of a matched antibody pair:

MP01008-1: 84086-2-PBS capture and 84086-4-PBS detection (validated in Cytometric bead array)

MP01008-2: 84086-2-PBS capture and 84086-3-PBS detection (validated in Cytometric bead array)

MP01008-3: 84086-2-PBS capture and 84086-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01008-1, CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84086-2-PBS. Detection antibody: 84086-4-PBS. Standard: Eg32087. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01008-2, CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84086-2-PBS. Detection antibody: 84086-3-PBS. Standard: Eg32087. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01008-3, CD244 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84086-2-PBS. Detection antibody: 84086-1-PBS. Standard: Eg32087. Range: 0.625-80 ng/mL