

For Research Use Only

DNA2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84105-3-PBS



Basic Information

Catalog Number:

84105-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001080449

GeneID (NCBI):

1763

UNIPROT ID:

P51530

Full Name:

DNA replication helicase 2 homolog (yeast)

Calculated MW:

120 kDa

Purification Method:

Protein A purification

CloneNo.:

241073B10

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Product Information

84105-3-PBS targets DNA2 as part of a matched antibody pair:

MP01042-1: 84105-3-PBS capture and 84105-1-PBS detection (validated in Cytometric bead array)

MP01042-2: 84105-3-PBS capture and 84105-2-PBS detection (validated in Cytometric bead array)

MP01042-3: 84105-4-PBS capture and 84105-3-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP01042-3, Human DNA2 Recombinant Matched Antibody Pair - PBS only. 84105-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY00221. 84105-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-2500 pg/mL



The mean DNA2 concentration was determined to be 2.00 ng/mL in Caco-2 cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP01042-1, DNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84105-3-PBS. Detection antibody: 84105-1-PBS. Standard: Ag17119. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01042-2, DNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84105-3-PBS. Detection antibody: 84105-2-PBS. Standard: Ag17119. Range: 0.781-100 ng/mL