

For Research Use Only

# SAMD9L Recombinant antibody, PBS Only (Capture)

Catalog Number: 84109-3-PBS



## Basic Information

<b>Catalog Number:</b> 84109-3-PBS	<b>GenBank Accession Number:</b> BC127118	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 219285	<b>CloneNo.:</b> 241250B1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8IVG5	
<b>Isotype:</b> IgG	<b>Full Name:</b> sterile alpha motif domain containing 9-like	
<b>Immunogen Catalog Number:</b> AG18541	<b>Calculated MW:</b> 1584 aa, 185 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

84109-3-PBS targets SAMD9L as part of a matched antibody pair:

MP01024-1: 84109-3-PBS capture and 84109-4-PBS detection (validated in Cytometric bead array)

MP01024-3: 84109-3-PBS capture and 84109-5-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01024-1, SAMD9L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84109-3-PBS. Detection antibody: 84109-4-PBS. Standard: Ag18541. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01024-3, SAMD9L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84109-3-PBS. Detection antibody: 84109-5-PBS. Standard: Ag18541. Range: 0.625-80 ng/mL