

For Research Use Only

Mouse CD80 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84137-2-PBS



Basic Information

Catalog Number: 84137-2-PBS	GenBank Accession Number: NM_009855.2	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 12519	CloneNo.: 241371C5
Source: Rabbit	UNIPROT ID: Q00609-1	
Isotype: IgG	Full Name: CD80 antigen	
	Calculated MW: 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Product Information

84137-2-PBS targets CD80 as part of a matched antibody pair:

MP01051-1: 84137-2-PBS capture and 84137-4-PBS detection (validated in Cytometric bead array)

MP01051-3: 84137-2-PBS capture and 84137-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01051-1, MOUSE CD80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84137-2-PBS. Detection antibody: 84137-4-PBS. Standard: Eg0840. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01051-3, MOUSE CD80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84137-2-PBS. Detection antibody: 84137-3-PBS. Standard: Eg0840. Range: 0.313-40 ng/mL.